


# Kistler MARS

Type 2875A...

## Measurement, Analysis and Reporting Software

Kistler MARS is innovative, comprehensive and user-friendly software for the complete analysis of force plate measurements. It supports routine diagnostics and research work in biomechanics, performance analysis, motor control behavior, rehabilitation medicine and other related fields.

- 20 different analysis modules
- User-friendly and intuitive
- Complete analysis of force measurements
- Uni- und bilateral Analysis
- Extensive data management
- Efficient preparation of reports



### Description

Kistler MARS is a routine diagnostics and research software for Kistler force plates. The 'Measurement, Analysis and Reporting Software' (MARS) supports Kistler data acquisition systems Types 5691A... and 5695B...

It analyzes the acquired force plate signals, calculates a large range of significant parameters, and provides graphical presentations of the measurements.

The software is based on a management unit that provides operational functions (add, edit, delete, assign, search, filter) to structure the data (projects, visits and subjects). The data is stored and managed in a database where it is accessible for comparison and reporting.

In each of the 20 different measurement modules the data is analysed for relevant output parameters. The input parameters and acquisition setup can be edited for each test separately.

The individual software functions are arranged intuitively and easy-to-use. All functions are well supported with extensive help information including how-to examples and with online video trainings on the MARS website.

### Application

Kistler MARS can be used for the evaluation of human movement such as static balance (Body Sway), dynamic balance (Tracking Shapes, Limits of Stability, etc.) locomotion and body transfer (Step Analysis, Forward Lunge, etc.), fast alternating movements (Tapping, Stamping, etc.) and strength and power (all vertical jumps, etc.). MARS additionally enables the user to carry out a free acquisition of signals coming out of the force plate, perform analysis and export data in different formats (raw signal, signal graph and parameters values).

The software calculates the standard parameters and many more evidence-based parameters to provide the details immediately. Detailed conclusions about muscular imbalances are possible when using two force plates.

Kistler MARS includes 20 modules to evaluate physical performance in the fields of power, anaerobic endurance, coordination and balance (see table page 2).

Name of the Test Module	Short Description	Full Version	Power Version	Balance Version
<b>Squat Jump</b>	Vertical jump test of concentric power for the lower extremities.	X	X	
<b>Counter Movement Jump</b>	Vertical jump test of eccentric-concentric power for the lower extremities.	X	X	
<b>Drop Jump</b>	Vertical drop jump test of eccentric-concentric power for lower legs. Testing is performed using progressively higher drop heights.	X	X	
<b>Jumps with Additional Weights</b>	3 consecutive vertical concentric and eccentric-concentric jumps. The test is performed using progressive loading with weights.	X	X	
<b>Repetitive Counter Movement Jumps</b>	Vertical jump test of endurance in eccentric-concentric conditions for lower extremities.	X	X	
<b>Repetitive Hopping</b>	Vertical jump test of endurance in eccentric-concentric conditions for lower legs.	X	X	
<b>Long Jump</b>	Situational horizontal jump test of eccentric-concentric power for lower extremities.	X	X	
<b>Squat</b>	Vertical movement test of concentric power for the lower extremities.	X	X	
<b>Stamping</b>	Test of maximal frequency and endurance of stamping for lower and upper extremities.	X	X	
<b>Tapping</b>	Test of maximal frequency, endurance and accuracy of tapping for lower and upper extremities.	X		X
<b>Forward Lunge</b>	Test for strength, good range of motion, balance and coordination.	X		X
<b>Sit-To-Stand</b>	Clinical test where the subject needs to rise from a seated to a standing position.	X		X
<b>Turn</b>	Clinical test, where the subject has to make two forward steps and then quickly turn for 180°.	X		X
<b>Step Analysis</b>	Situational test of vertical, anterior-posterior and medio-lateral load of the lower extremity during locomotion.	X		X
<b>Body Sway</b>	Test of body sway during sustaining static posture (quiet stance or any other). Typical applications: sport performance, lower extremity injuries, brain concussion, etc.	X		X
<b>Tracking Shapes</b>	A set of tests which involve precise active tracking of the centre of pressure way defined by a matrix shape that is displayed on a screen as a feed back to the subject.	X		X
<b>Tracking Curves</b>	A set of tests which involve precise active tracking of the centre of pressure way defined by a matrix curve that is displayed on a screen as a feed back to the subject.	X		X
<b>Limits of Stability</b>	A test of maximal range of voluntary body leaning in different directions.	X		X
<b>Landing</b>	Dynamic balance test of postural stability for different types of landing.	X		X
<b>Symmetry</b>	Test of postural symmetry in weight bearing function of the lower extremities during upright stance and (semi)squat.	X		X

## Application Examples and Screenshots




Fig. 1: Main window with management tools on the left side and module selection on the right



Fig. 2: The Balance Version provides all modules related to locomotion, movement control and balance




Fig. 3: Example: Counter Movement Jump with automatic event detection




Fig. 4: Example: Squat Jump result table




Fig. 5: Example: Shape Tracking

**Included Accessories**

- License on USB Key

**Optional Accessories**

- Kistler force plates
- Kistler DAQ systems

Type/Art. No.  
9260...  
9286...  
9281...  
9287...  
5691A...  
5695B...

**Ordering Key**

Type 2875A ☐

**Kistler MARS**

Full Version Measurement, Analysis and Reporting Software	1
Power Software Measurement, Analysis and Reporting Software Power modules only	3
Balance Software Measurement, Analysis and Reporting Software Balance modules only	4

The software is developed by S2P.

Copyright by **S2P**  
SCIENCE TOP PRACTICE www.s2p.si

This information corresponds to the current state of knowledge. Kistler reserves the right to make technical changes. Liability for consequential damage resulting from the use of Kistler products is excluded.

©2012 ... 2014, Kistler Group, Eulachstrasse 22, 8408 Winterthur, Switzerland  
Tel. +41 52 224 11 11, Fax +41 52 224 14 14, info@kistler.com, www.kistler.com  
Kistler is a registered trademark of Kistler Holding AG.