

# Multicomponent Force Plate

Type 9281E...

## for Dynamic Applications in Biomechanics, $F_z$ -10 ... 20 kN

Multicomponent force plate with wide range for measuring ground reaction forces, moments and the center of pressure in biomechanics.

- Extremely wide measuring range
- Excellent measuring accuracy
- High natural frequency
- Versatile
- Threshold  $F_z < 250$  mN

### Description


The multicomponent force plate Type 9281E... consists of a 600x400 mm aluminum sandwich top plate of advanced, lightweight construction and four built-in piezoelectric 3-component force sensors. Thus it is extremely rigid overall, and allows measurements over a very wide useful frequency range.

Thanks to the special properties of the piezoelectric sensors, the plate is highly sensitive and can simultaneously measure very dynamic phenomena involved in a wide range of applications.

### Application

This force plate is designed specifically for use in basic research and sport. Its extensive range and high rigidity allow it to be employed across a very wide spectrum of measuring tasks and application sectors. Despite the very generous measuring range of -10 ... 20 kN, it offers excellent accuracy and linearity, and even under a large preload allows precise measurement of minute forces. In all these situations the plate can be mounted in any position without affecting the measurement result in any way.

The Type 9281EA has a built-in charge amplifier compatible with all of the common motion analysis systems.



### Technical Data

Dimensions		mm	600x400x100
Measuring range	$F_x, F_y$	kN	-10 ... 10
	$F_z$	kN	-10 ... 20
Overload	$F_x, F_y$	kN	-15/15
	$F_z$	kN	-10/25
Linearity	% FSO		$< \pm 0,2$
Hysteresis	% FSO		$< 0,3$
Crosstalk	$F_x \leftrightarrow F_y$	%	$< \pm 1,5$
	$F_x, F_y \rightarrow F_z$	%	$< \pm 1,5$
	$F_z \rightarrow F_x, F_y$	%	$< \pm 0,5^{1)}$
Rigidity	x-axis ( $a_y = 0$ )	N/ $\mu$ m	$\approx 250$
	y-axis ( $a_x = 0$ )	N/ $\mu$ m	$\approx 400$
	z-axis		
	( $a_x = a_y = 0$ )	N/ $\mu$ m	$\approx 30$
Natural frequency	$f_n$ (x, y)	Hz	$\approx 1\,000$
	$f_n$ (z)	Hz	$\approx 1\,000$
Operating temperature range		$^{\circ}$ C	0 ... 60
Weight		kg	16
Degree of protection	EN 60529:1992		IP65

<sup>1)</sup> inside sensor rectangle

**Force Plate with Built-in 8-Channel Charge Amplifier, Type 9281EA**

Calibrated range	$F_x, F_y$	kN	0 ... 5
	$F_z$	kN	0 ... 20
Calibrated partial range	$F_x, F_y$	kN	0 ... 1,25
	$F_z$	kN	0 ... 5
Sensitivity range 1	$F_x, F_y$	mV/N	$\approx 40^{2)}$
	$F_z$	mV/N	$\approx 18^{2)}$
Sensitivity range 4	$F_x, F_y$	mV/N	$\approx 2,0^{2)}$
	$F_z$	mV/N	$\approx 0,9^{2)}$
Ratio ranges 1:2:3:4			1:5:10:20 <sup>3)</sup>
Threshold		mN	$<250^4)$
Drift		mN/s	$<\pm 10$
Supply voltage		VDC	10 ... 30
Supply current		mA	$\approx 45$
Output voltage		V	0 ... $\pm 5$
Output current		mA	-2 ... 2
Control inputs (optocoupler)		V	5 ... 45
		mA	0,4 ... 4,4

**Force Plate without Charge Amplifier, Type 9281E**

Calibrated range	$F_x, F_y$	kN	0 ... 10
	$F_z$	kN	0 ... 20
Calibrated partial range	$F_x, F_y$	kN	0 ... 1
	$F_z$	kN	0 ... 2
Threshold	$F_x, F_y, F_z$	mN	$<50$
Sensitivity	$F_x, F_y$	pC/N	$-7,5^{2)}$
	$F_z$	pC/N	$-3,8^{2)}$

- <sup>2)</sup> nominal value
- <sup>3)</sup>  $\pm 0,5$  % accuracy
- <sup>4)</sup> only range 1

Conforms to the **CE** safety standards (73/23/EG) for electrical equipment and systems:

EN 60601-1:2005, EN 61010-1:2001

and the EMC standards (89/336/EG):

EN 60601-1:2005 (EN 55022 Class B), EN 61000-6-3:2004

(EN 55022 Class B), EN 61000-6-4:2001 (EN 55011 Class B),

EN 60601-1:2005, EN 61000-6-1:2001, EN 61000-6-2:2005

**Dimensions**




Fig. 1: Dimensions of multicomponent force plate Type 9281E...

9281E\_000-711e-02\_14

**BioWare®**

BioWare software is the engine behind the force plate system. It collects data from the force plates, converts the trials into useful information and plots the results. The force plates and charge amplifiers are fully remote controlled by BioWare thus making the system extremely flexible and easy-to-use.

**Parameters of Gait**




Fig. 2: Ground reaction forces (GRF)




Fig. 3: Center of pressure (COP)




Fig. 4: Frictional torque  $T_z$




Fig. 5: Force vector

**Other functions**

- Coefficient of friction (COF)
- Frequency analysis, statistics, digital filters
- Full Windows® functionality

Windows® is a registered trade mark of Microsoft Corporation.

This information corresponds to the current state of knowledge. Kistler reserves the right to make technical changes. Liability for consequential damage resulting from the use of Kistler products is excluded.

BioWare provides several performance specific evaluations.

**Parameters of Countermovement Jump CMJ**




Fig. 6: Jump force




Fig. 7: Power




Fig. 8: Jump height (COM)




Fig. 9: Force gradient (Explosivity)

**Other parameters**

- Acceleration, velocity and displacement of the center of mass (COM)
- Work, energy, impulse
- Statistics, digital filters

9281E\_000-711e-02.14

**Typical Measuring Chains**





			
Force plate with charge amplifier Type 9281EA	Connection cable Type 1759A...	DAQ system (USB 2.0) Type 5691A1	Laptop (provided by user) with BioWare software

Fig. 10: Configuration of a typical measuring chain with Kistler DAQ system BioWare®




			Ch 1 = $F_x$ 1+2    Ch 5 = $F_z$ 1 Ch 2 = $F_x$ 3+4    Ch 6 = $F_z$ 2 Ch 3 = $F_y$ 1+4    Ch 7 = $F_z$ 3 Ch 4 = $F_y$ 2+3    Ch 8 = $F_z$ 4
Force plate with charge amplifier Type 9281EA	Connection cable Type 1757A...	External Control Unit (8xBNC neg.) Type 5233A2	DAQ system provided by user (8 analog channels)

Fig. 11: Configuration of a typical measuring chain with DAQ system provided by user

**Included Accessories**

**For Type 9281E...**

- |  |                        |
|--|------------------------|
| • 1 Set shims                                | 7.050.011              |
| • 4 Eye bolts M6<br>with washer              | 6.170.007<br>6.220.040 |
| • 4 Hexagon socket head cap screws<br>M12x25 | 6.120.106              |
| • 1 Hexagon socket wrench                    | 1391                   |
| • 1 Voltage equalizing cable                 | 5.590.175              |
| • 4 Installation handles                     | 7.511.437              |

**Optional Accessories**

**For Type 9281EA with built-in charge amplifier**

- |  |              |
|--|--------------|
| • 16ch DAQ-System for BioWare (USB 2.0)  | 5691A1       |
| • Connection cable for 5691A, angle plug | 1759A...     |
| • 64ch DAQ-System for BioWare (USB 2.0)  | 5695B1       |
| • Connection cable for 5695B, angle plug | 1700A105A... |
| • External control unit (BNC out)        | 5233A2       |
| • Connection cable for Type 5233A...     | 1757A...     |
| • DAQ system BioWare (PCI-Bus)           | 2812A...     |

**For Type 9281E with charge output**

- |  |          |
|--|----------|
| • External charge amplifier              | 9865E... |
| • Connection cable, angle plug connector | 1686A... |
| • DAQ system BioWare (PCI-Bus)           | 2812A... |

**Mounting frame for Type 9281E...**

- |   |            |
|---|------------|
| • Standard mounting frame                             | 9423       |
| • Other mounting frames for multiple<br>installations | on request |

**Ordering Key**

**Multicomponent Force Plate**

with charge output	-
with built-in charge amplifier	A

Type 9281E



9281E\_000-711e-02\_14

BioWare® is a registered trade mark of Kistler Holding AG.